



MMSlog® MODULAR AUXILIARY TROLLEY

OPENED MODULAR AUXILIARY TROLLEY WITH 4, 5 & 8 GUIDE RAILS FOR ISO3394 TRAYS & BASKETS

- Opened auxiliary trolley with standard dimensions 600*400 mm
- Structure made of epoxy-coated steel profile
- Trolley body in square profile 20*20 mm
- Tubular lateral push handle 20 mm diameter
- Built-in brackets on top to support the HPL working surface, dimensions 600*400*6 mm
- Intern span fitted with 4, 5 or 8 pairs of lateral guides rails, for ISO3394 trays and/ or baskets
- Trolley base equipped with 4 rotary wheels diameter 75 mm, 2 of them with brake

Reference	External dimensions			Number of guides	Pairs	Wheels
	Product Code	Height	Width	Depth	Body	Guides
CA105000000	825 mm	755 mm	475 mm	1	5	75 mm Diam. 4 Rotary Wheels 2 of them with brake
CA106000000	945 mm	755 mm	475 mm	1	6	
CA210000000	825 mm	995 mm	675 mm	2	5+5	
CA212000000	945 mm	995 mm	675 mm	2	6+6	